PTS Panels® Multi-Chemistry Controls

for use with CardioChek® Brand Analyzers. Not suitable for self-testing.
For Professional Use Only.

PS-000597EFC Rev. 1 06/11

INTENDED USE

Use these controls with PTS Panels Test Strips as quality control materials to check for proper performance of your analyzer. These controls are intended for use in estimating precision and detecting systematic analytical deviations that may arise from reagent or analytical instrument variation. They are intended for use by healthcare professionals in both physicians' offices and in acute and convalescent care facility bedside testing. Currently PTS, Inc. tests for cholesterol and triglycerides are licensed in Canada. This product is only intended for use with those tests in Canada.

SUMMARY AND EXPLANATION

PTS Panels Multi-Chemistry Controls contain cholesterol, triglycerides, glucose, and ketone, which react with test strips to produce color. The darker the color, the higher the analyte concentration. The result for each level of control should be in the range specified on the Quality Control Range Card included with these materials. Controls should be run to verify the test system performance, when results are questionable or to comply with your facility's quality control requirements. (Reference analyzer User Guide)

CHEMICAL COMPOSITION

PTS Panels Multi-Chemistry Controls contain glucose, glycerol, cholesterol, DL-β-Hydroxybutyric Acid, MIT (methyl isothiazolidone) and inert ingredients.

PRECAUTIONS

- Do not eat or drink contents.
- For *in vitro* diagnostic use.
- Avoid contact with eyes or skin.

STORAGE AND HANDLING

Product is stable at room temperature until marked expiration date. Do not freeze. Optimal long term storage is 2-8°C. Product may be stored for up to four months at room temperature (20-30°C) once opened. Product is ready to use and is clear. Discard if appears cloudy. Recap after using. Use care not to contaminate dropper tip. Do not use beyond the expiration date on the vial.

PROCEDURE

Materials provided:

Two (2) vials of Multi-Chemistry Control Materials (one each of Level 1 and Level 2)

Materials required, but not provided:

- CardioChek Brand Analyzer
- PTS Panels Test Strips and MEMO Chip® of same lot number

PLEASE READ INSTRUCTIONS CAREFULLY BEFORE TESTING

Test Procedure:

Please refer to analyzer User Guide for the proper procedure on how to run a quality control test. The following procedure applies to the CardioChek brand analyzers:

- Insert MEMO Chip that matches the lot of strips.
- Press one of the two buttons to turn the instrument ON.
- Press the NEXT button until the display reads UTILITY.
- Press the ENTER button.
- Press the NEXT button until RUN CONTROL is displayed.
- At the RUN CONTROL analyzer display option, press ENTER.
- Insert the test strip into analyzer.
- Wait until analyzer displays APPLY SAMPLE.
- Invert to mix control.
- Remove cap and turn bottle upside down straight above test strip sample application window.
- Carefully, squeeze vial to deliver one drop of control solution on the test strip.
- Replace control vial cap.
- Results will be displayed in about 2 minutes.

EXPECTED RESULTS

Results for the control should be in the range specified on the Quality Control Range Card. The control results are also available on the PTS, Inc. website: www.cardiochek.com. If test results are above or below the range specified, check to see if the test strips or Controls are expired. If they are expired, test again with unexpired strips and controls. Clean the analyzer before retesting (See User Guide). For customer service and/or technical support call PTS, Inc. Customer Service at +1-877-870-5610 (Toll-free, inside US) or +1-317-870-5610.
email: inforequest@cardiochek.com

Authorized European Representative
per IVDD 98/79/EC
MDSS Schiffgraben 41
30175 Hannover, Germany

Polymer Technology Systems, Inc.
Indianapolis, Indiana 46268 USA

© Copyright 2011 Polymer Technology Systems, Inc.

Use By

REF Catalog number

Batch Code

Manufacturer

Temperature limitation

In vitro diagnostic medical device

Consult instructions for use

This product fulfills the requirements of Directive 98/79 EC for *in vitro* diagnostic medical devices.

EC REP Authorized representative in the European Community

CONTROL Control

PTS Panels® Multi-Chemistry Controls

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