

**RUN TEST**

( [arrow] to select)

USE CODE #

INSERT STRIP



**MEMORY**

( [arrow] to select)

[arrow] CHEM  
( [arrow] to select)

- [arrow] CHOL
- [arrow] GLU
- [arrow] HDL
- [arrow] TRIG
- [arrow] LDL
- [arrow] KET
- [arrow] LIPIDS
- [arrow] EXIT

[arrow] CONTROL

[arrow] EXIT - - - [arrow] return to MEMORY

[arrow] return to CHEM



**UTILITY**

( [arrow] to select)

[arrow] CK STRIP

[arrow] RUN CONTROL

[arrow] CLEAR MEMORY

[arrow] EXIT - - - [arrow] return to UTILITY

[arrow] return to CK STRIP



**USER SET**

( [arrow] to select)

[arrow] PRINT

[arrow] SOUND

[arrow] SET TIME

[arrow] SET DATE

[arrow] UNITS\*

[arrow] EXIT - - - [arrow] return to USER SET

[arrow] return to PRINT

\* Units may be locked, and therefore UNITS will not appear.